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ENGLISH

INDICATIONS

For the support of restorations and prosthetic crowns.

COMPOSITION

Peeso drills.....Stainless steel
Rubber Stops.....Latex
Template for radiographic measurement.....Plastic

REFORPOST® GLASS FIBER: Glass fiber (80%); pigmented resin (19%), stainless steel filament (1%).
REFORPOST® CARBON FIBER: Carbon fiber (80%); pigmented resin (19%), stainless steel filament (1%).

TECHNIQUE OF USE

A. ROOT CANAL PREPARATION

- Select the pin using the template and the periapical x-ray of the tooth (photo 1);
- Remove some of the root canal filling with heated instruments to the predetermined depth. The use of drills can cause apical displacement of the gutta-percha. Allow at least 3 to 5 mm of filling at the root apex due to the presence of accessory canals, which must remain sealed to prevent bacterial contamination (photo 2);
- Prepare the canal with the drill corresponding to the diameter of the selected post (photo 3).

B. PREPARATION OF REFORPOST FOR CEMENTATION

- Position the pin to verify its adaptation. Make a mark on the pin 2 to 3 mm from the antagonist tooth (photo 4) – and cut it with a diamond bit at high rotation under cooling. Do not use scissors or pliers, because the maceration of the fibers changes the mechanical properties of the pin;
- Clean the post with alcohol and dry it (photo 5);
- Apply a layer of silane (SILANO ANGELUS®) and let it dry for 1 minute (photo 6);

Ⓛ WARNING: It is not necessary to apply the silane on REFORPOST® CARBON FIBER.

4. Apply the ANGELUS® FUSION-DURALINK® CATALYST to the post according to its instructions (photo 7).

C. CEMENTATION

- Clean the canal with the air-water spray and dry it with air and absorbent paper points (photo 8);

- Etch the canal with ANGELUS® 37% PHOSPHORIC ACID for 15 seconds (photo 9);
- Wash again the canal and dry it with air and absorbent paper points (photo 10). Apply ANGELUS® FUSION-DURALINK® PRIMER in the canal and wait for 1 minute;
- Apply ANGELUS® FUSION-DURALINK® CATALYST according to the instructions. The adhesive components must not be mixed before application. They are polymerized chemically without light curing (photo 11);
- Spatulate the resinous cement selected* and fill the conduit;
- Put the post in the canal, remove excess of the cement and wait for the polymerization. If the cement is dual, light cure it (photo 13).

Ⓛ Zinc phosphate and glass ionomer cements (conventional and resin-modified) can also be used, but their mechanical properties are inferior to those of the resin cements. Auto-adhesive cements eliminate the needs for the steps of acid conditioning and adhesive application.

D. CROWN BUILD-UP

- Do another etching with ANGELUS® 37% PHOSPHORIC ACID and apply an adhesive on the remaining tooth structures and post. Build the crown core with composite resin (photos 14 and 15), or with ANGELUS® REFORCORE®. REFORCORE® is a prefabricated glass fiber core, for the technique of direct cores (photos 16 and 17).

E. USE OF ACCESSORY POSTS IN WEAKENED ROOTS AND WIDE CANALS

- Fill the remaining spaces of the canal with ANGELUS® REFORPIN® accessory posts. The decrease in thickness of the cement and the better retention provided by REFORPIN® reduce the risk of post/core loosening (photos 18 and 19).

WARNINGS | PRECAUTIONS

- Use of direct cores with fiber posts is not indicated for extensive oral rehabilitation in patients with parafunctional habits;
- Use of the direct core with fiber post should only be indicated to teeth with a minimum height of 2 mm of coronal dentin;
- Fiber posts should only be cut with water refrigerated high speed diamond burs. Do not use scissors or cutting pliers. Use goggles and a mask to prevent inhalation of particles generated during cutting;
- Fiber posts must be completely covered by the cementing agent or the composite resin. Exposed fibers absorb water which deteriorates the mechanical properties of the posts.

Ⓛ IMPORTANT: The indication of any post depends on careful evaluation of the quality of endodontic treatment and periapical health.

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